



LED PX3 Device Quick User Guide



Product Name	LED PX3 Device
Catalog #	AT8003-D
Outside Dimensions	7 inch x 5 1/2 inch x 6 1/4 inch (L x W x H) 18 cm x 14 cm x 16 cm (L x W x H)
Device Operation	1. Power-On: The screen activates automatically upon powering on the device and enters time-setting mode by default. A session duration must be configured before the light can be activated. 2. Light Control: Press the button once to turn the light on. Press again to pause. Press and hold for more than 2 seconds to reset the timer to the currently configured duration. There is an automatic safety switch that turns off the UV when the device is opened. 3. Function Selector (Gear-Shaped Button): Press to cycle between three display modes: Time Setting, Light Intensity, and UV Testing. Use the + and – buttons to adjust the selected parameter – time adjusts in 5-minute increments, intensity in 5% increments, and UV units between $\mu\text{W}/\text{cm}^2$ and mW/cm^2. Note: UV output measurement requires the optional UV sensor accessory, available for purchase separately. Without the accessory connected, the UV Testing screen will not display a value. 4. Cooling Fan: The fan activates automatically when the light turns on and continues to run for 10 seconds after the light turns off.

For more information, please visit www.alphathera.com or contact support@alphathera.com.

AlphaThera, Inc.

Version 1.1 September 2026

Scan to watch the
How To Video:

