



oYo-Link® LED PX2 Device

Product Specification Sheet

Product Name	oYo-Link® LED PX2 Device
Catalog #	AT8001
User Manual:	See the LED PX2 Device user manual at: https://alphathera.com/user-manuals
Outside Dimensions	8 / inch x 9 / inch x 4 inch (L x W x H) 22cm x 27cm x 10cm
Voltage	Input: 100-240V AC 50/60Hz Output: DC 14V/2A
Fuse	3A
Bulb Type	LED bulb 365 nm (9W x2 standard)
Bulb Base	Bulb Base G23
Power	18W
Timer	60 mins and 120 mins; Use 120 mins (2 hours) to ensure maximum conjugation efficiency
Fan	Yes
Notes	<p>This LED device is only to be used for in-vitro laboratory research and not to be used for human or clinical samples. Avoid working in moist or humid environments. Always disconnect the main power cord before cleaning the unit or replacing the fuse or LED bulbs. WARNING: The blacklight bulbs used in this device emit long wavelength ultraviolet (UV) light (365 nm), which can pose potential health risks including injuries to the skin and eyes. Please exercise caution when the device is running and avoid eye or skin exposure.</p>

For more information, please visit www.alphathera.com or contact support@alphathera.com.